



BAR AND CLIP BOXES AND SIZING LINERS

The must have accessories for busy clubfoot clinics needing to operate at maximum efficiency.

Bar and Clip Boxes organise a complete stock of ADM Modular System Bar and Clip components making it easy to assemble a fully configured External Rotation Bar to your patient's requirements in just a few seconds. Sizing Liner sets are the best and fastest way of determining which size orthosis is required for your patient.



Clubfoot brace fitting made easy:

- Measure patient's shoulder width
- Take required bar and clips from box and click together
- Use Sizing Liners to determine the best AFO size for each foot
- Fit configured brace to the patient

No adjustments, no screws, no allen keys, no fuss

Small Bar & Clip Boxes



- 16 clips (8 left & 8 right)
- 3 x complete sets of bars (21 bars)
- 6 x clip removal tools
- Measuring tape
- Instructions for Use
- Contents card

Large Bar & Clip Boxes



- 8 clips (4 left & 4 right)
- 3 x complete sets of bars (18 bars)
- 6 x clip removal tools
- Instructions for Use
- Contents card

Sizing Liners



- Small: 9 liner sizes left and right (0000 to 5)
- Large: 7 liner sizes (6 to 12)
- Liners organised with removal cord and identification tags

Code	Product Description
BCBS70604530	ADM Bars and Clips Box Small 70, 60, 45, 30 DEG
BCBS604530	ADM Bars and Clips Box Small 60, 45, 30 DEG
BCBS704530	ADM Bars and Clips Box Small 70, 45, 30 DEG

Note: BCBS70604530 has 2 pairs of each clip type. BCBS604530 has 4 pairs of 60 degree clips, but no 70 degree clips. BCBS704530 has 4 pairs of 70 degree clips, but no 60 degree clips.

Ordering replacement components: Download our our easy order form from www.c-prodirect.com

Code	Product Description
BCBL70604530	ADM Bars and Clips Box Large 70, 60, 45, 30 DEG
BCBL604530	ADM Bars and Clips Box Large 60, 45, 30 DEG
BCBL704530	ADM Bars and Clips Box Large 70, 45, 30 DEG

Note: BCBL70604530 has 1 pair of each clip type. BCBL604530 has 2 pairs of 60 degree clips, but no 70 degree clips. BCBL704530 has 2 pairs of 70 degree clips, but no 60 degree clips.

Ordering replacement components: Download our our easy order form from www.c-prodirect.com

Code	Product Description
SLNR	Small Liner Set sizes 0000 to 5 Left and Right
LLNR	Large Liner Set sizes 6 to 12 Left and Right

Note: Sizing Liners are identical to the liners of ADM AFOs and are supplied labelled and joined into sets with a cord. A heel port enables the position of the heel to be seen.

Price: Please contact C-Pro Direct or your distributor for price information.